



Accessories

Side Rail & Bed End Bumpers

To avoid entrapment Side Rail Bumpers are available to fit different side rails. The bumpers are fastened using Velcro, ensuring that the bed environment is safe.

The bumpers are made from a padded, vapour permeable wipe clean fabric which is machine washable up to 60°C.

Bed End Bumpers are available to the same specification if required.



Wedge Set – Side Wedge set with draw sheet

This method of preventing a patient from falling or rolling out of bed requires no use of side rails. It comprises of a draw sheet which easily attaches to the mattress with quick release buckles and two wedges which are positioned along each side of the mattress using Velcro fastenings.

The draw sheet and the wedges have covers which are made from a vapour permeable wipe clean fabric. The draw sheet can also be machine washed up to 60°C if required. A cotton draw sheet is also available on request.

Wedge size 90cm long x 20cm high, bespoke sizes available on request.



Tumble Mats

Tumble Mats are available to help ensure safety around beds, both standard and low level.

Made from high density foam and covered in a non-slip two way stretch vapour permeable, wipe clean fabric. The covers can be machine washed up to 60°C if required.

The mats fold away easily for storage.

Standard size: 198 x 60 x 5 cm

Double thickness and bespoke sizes available on request.

Tumble Mat alongside bed



Tumble Mat folded for storage



SCAN MOBILITY LTD

Marsh House, Ralph's Wife's Lane,
Banks, Lancashire, PR9 8ER

Tel: **01704 233344** Fax: 01704 233355

E-mail: sales@scanmobility.co.uk www.scanmobility.co.uk