


SERENITY HYBRID

1 in 3 Active Mattress from SCAN MOBILITY

Hybrid support surface combine foam and air to maximize the benefits offered by both static and alternating surfaces. This simplifies choice, as one piece of equipment is suitable across a much broader range of patients. This system can be powered and non-powered giving a much wider range of pressure relief with one mattress.

Cover Technology The mattress cover is made from a new revolutionary versatile fabric, which is multi stretch and breathable. This helps to contribute to the reduction of pressure sore development by helping to prevent a build-up of body moisture, but at the same time remaining fluid proof. This therefore reduces heat build-up, meeting stringent infection control requirements.

The cover offers improved durability and has 10,000-ppm chlorine resistance as well as conforming to BS7175.



OEKI-TEX 100 compliant



The Serenity Hybrid comprises a high-specification foam mattress made from castellated foam. The foam is encased in air cells, which can be removed to make drop zones.

Should the condition of the patient deteriorate, the pump can be used to activate the air cells around the foam for alternating pressure relief, therefore reducing pressure ulcer incidence and cost.

The Pump is a digital powered unit, providing the mattress with a **1 in 3** cycle giving greater patient support.

The **1 in 3** cycle provides excellent support whilst achieving high-pressure relief with quiet operation.

Lock out mode is also included.

The unit works on **low-pressure therapy**. This allows greater immersion into the mattress, helping to provide optimum pressure distribution making it good for very high-risk patients.



Dimensions

Short 183x90x20cm

Standard 200x90x20cm

Max Patient Weight

39st/250kg

Control Unit Weight

2.2kg

Control Unit Dimensions

W29cmxH16.5cmxD11cm

Risk Category

Very High

2-YEAR WARRANTY